



Cartridges



MODEL	CAR-20	CAR-20H	CAR-30	CAR-40	CAR-50	CAR-60	
Type	Moving Coil	High output MC	Moving Coil	Moving Coil	Moving Coil	Moving Coil	
Coil Wire	STD Copper	STD Copper	5N Copper	4N Silver	4N Silver	4N Silver	
Cantilever Material	Aluminium	Aluminium	Boron	Boron	Sapphire	Diamond	
Stylus	Synthetic Elliptical	Synthetic Elliptical	Microridge	Microridge	Microridge	Microridge	
Typical figures	Frequency Response	10 Hz - 33 kHz	10 Hz - 33 kHz	10 Hz - 35 kHz	10 Hz - 40 kHz	10 Hz - 45 kHz	10 Hz - 45 kHz
	Output Voltage (3.54 cm/1 kHz)	0.3 mV	2 mV	0.3 mV	0.3 mV	0.3 mV	0.3 mV
	Channel Balance (1 kHz)	<1 dB	<1.5 dB	<1 dB	<1 dB	<0.5 dB	<0.5 dB
	Channel Separation (1 kHz)	>23 dB	>21 dB	>25 dB	>28 dB	>30 dB	>30 dB
	Tracking Force	2.0 gr	2.0 gr	2.0 gr	2.0 gr	2.0 gr	2.0 gr
	Compliance	8×10^{-6} cm/dyne	8×10^{-6} cm/dyne	10×10^{-6} cm/dyne	10×10^{-6} cm/dyne	10×10^{-6} cm/dyne	10×10^{-6} cm/dyne
	Trackability	>70 μ m/2.0 gr	>60 μ m/2.0 gr	>70 μ m/2.0 gr	>70 μ m/2.0 gr	>70 μ m/2.0 gr	>70 μ m/2.0 gr
	Internal Impedance	4 Ω	130 Ω	4 Ω	6 Ω	6 Ω	6 Ω
Load Impedance	<100 Ω	47 k Ω	<100 Ω	<100 Ω	<100 Ω	<100 Ω	
Net Weight	17 gr	17 gr	17 gr	17 gr	17 gr	17 gr	



Accessories



RD Ultrasonic kit

RD kit will allow you to submerge records in the ultrasonic cleaning bath and position them in a suitable drying position.



Ebony clamp



Outer clamp



Headshell container single



Headshell container 4

